



Product: Fmoc-Aib-OH  
SDS No: 231109  
Revision: 1.0

## Safety Data Sheet

### Section 1 - Product and Company Information

Product Name: Fmoc-Aib-OH  
Product Code: 231109  
Company: CSBio(Shanghai) Ltd  
EMERGENCY CONTACTS: Dr. Simon Zhu +86-21-50781655

### Section 2 - Hazards Identification

#### 2.1 Classification of the substance or mixture

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Acute toxicity, oral(Category 4) H302

Skin corrosion/irritation (Category 2) H315

Serious eye damage/eye irritation(Category 2A ) H319

Specific target organ toxicity, single exposure;Respiratory system(Category 3) H335

No Resource File

#### 2.2 GHS Label elements, including precautionary statements

Pictogram



Signal word

Warning

Hazard statement(s)

H302 Harmful if swallowed

H315 Causes skin irritation

H319 Causes serious eye irritation

H335 May cause respiratory irritation

Precautionary statement(s)

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

#### 2.3 Other hazards

Additional precautionary phrases are located throughout the safety data sheet.

### Section 3 - Composition/Information on Ingredients

Substance Name: Fmoc-Aib-OH  
CAS No. 94744-50-0  
Formula  $C_{19}H_{19}NO_4$

### Section 4 - First Aid Measures

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Aib-OH

SDS No: 231109

Revision: 1.0

#### 4.1 Description of first aid measures

##### General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area

##### If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

##### In case of skin contact

Take off contaminated clothing and shoes immediately. Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

##### In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician. Continue rinsing eyes during transport to hospital.

##### If swallowed

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

#### 4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

#### 4.3 Indication of any immediate medical attention and special treatment needed

No data available

### Section 5 - Fire Fighting Measures

#### 5.1 Extinguishing media

Suitable extinguishing media: alcohol-resistant foam, dry chemical or carbon dioxide

Unsuitable extinguishing media: no data

#### 5.2 Special hazards arising from the substance or mixture

In combustion toxic fumes may form.

#### 5.3 Advice for firefighters

Wear self-contained breathing apparatus.

Wear protective clothing to prevent contact with skin and eyes.

### Section 6 - Accidental Release Measures

#### 6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. Do not attempt to take action without suitable protective clothing.

For personal protection see section 8.

#### 6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains or river. Alert the neighbourhood to the presence of fumes or gas.

#### 6.3 Methods and materials for containment and cleaning up

Mix with sand or vermiculite. Sweep up and shovel. Transfer to a closable, labelled salvage container for

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Aib-OH

SDS No: 231109

Revision: 1.0

disposal by an appropriate method.

#### 6.4 Reference to other sections

For disposal see section 13.

### Section 7 - Handling and Storage

#### 7.1 Precautions for safe handling

Keep container tightly closed in a dry and well-ventilated area. Containers which are opened must be carefully resealed and kept upright to prevent leakage. Inert atmosphere, ( $\leq 8^{\circ}\text{C}$ )

#### 7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Cool dry place ( $\leq 8^{\circ}\text{C}$ )

Storage class (TRGS 510): Non Combustible Solids

#### 7.3 Specific end use(s)

The end use(s) have not been fully determined. The substance is supplied for Research and Development purposes by professionals only.

### Section 8 - Exposure Controls / PPE

#### 8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

#### 8.2 Exposure controls

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Control of environmental exposure

No special environmental precautions required.

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Aib-OH

SDS No: 231109

Revision: 1.0

## Section 9 - Physical/Chemical Properties

Appearance:	White crystalline powder
Molecular Weight:	325.4
pH:	N/A
BP/BP Range:	N/A
MP/MP Range:	182.0-188.0°C
Freezing Point:	N/A
Vapor Pressure:	N/A
Vapor Density:	N/A
Saturated Vapor Conc.	N/A
SG/Density:	N/A
Bulk Density:	N/A
Odor Threshold:	N/A
Volatile%:	N/A
VOC Content:	N/A
Water Content:	≤0.50%
Solvent Content:	N/A
Evaporation Rate:	N/A
Viscosity:	N/A
Surface Tension:	N/A
Partition Coefficient:	N/A
Decomposition Temp.:	N/A
Flash Point:	N/A
Explosion Limits:	N/A
Flammability:	N/A
Autoignition Temp:	N/A
Refractive Index:	N/A
Optical Rotation:	N/A
Miscellaneous Data:	N/A
Solubility:	N/A

## Section 10 - Stability and Reactivity

### 10.1 Reactivity

No data available

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Aib-OH  
SDS No: 231109  
Revision: 1.0

- 10.2 Chemical stability  
Stable under recommended storage conditions.
- 10.3 Possibility of hazardous reactions  
No data available
- 10.4 Conditions to avoid  
No data available
- 10.5 Incompatible materials  
Strong oxidizing agents
- 10.6 Hazardous decomposition products  
Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)  
Other decomposition products - No data available  
In the event of fire: see section 5

## Section 11 - Toxicological Information

- 11.1 Information on toxicological effects
- Acute toxicity  
No data available
- Inhalation: No data available
- Dermal: No data available
- No data available
- Skin corrosion/irritation  
No data available
- Serious eye damage/eye irritation  
No data available
- Respiratory or skin sensitisation  
No data available
- Germ cell mutagenicity  
No data available
- Carcinogenicity  
IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  
ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.  
NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.  
OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
- Reproductive toxicity  
No data available
- Specific target organ toxicity - single exposure  
No data available

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Aib-OH  
SDS No: 231109  
Revision: 1.0

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Additional Information

RTECS: Not available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

## Section 12 - Ecological Information

### 12.1 Toxicity

No data available

### 12.2 Persistence and degradability

No data available

### 12.3 Bioaccumulative potential

No data available

### 12.4 Mobility in soil

No data available

### 12.5 Results of PBT and vPvB assessment

No data available

### 12.6 Other adverse effects

No data available

## Section 13 - Disposal Considerations

### 13.1 Waste treatment methods

Disposal Operations

Consult state, local or national regulations for proper disposal. Hand over to authorised disposal company as hazardous waste.

Disposal of Packaging

Disposal must be made according to official regulations.

## Section 14 - Transport Information

DOT (US)

Classed non-hazardous for shipping.

IMDG

Classed non-hazardous for shipping.

IATA

Classed non-hazardous for shipping.

## Section 15 - Regulatory Information

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Aib-OH  
SDS No: 231109  
Revision: 1.0

Safety, health and environmental and national regulations:

Product is not subject to any additional regulations or provisions

Safety Assessment

No Chemical Safety Assessment

## Section 16 - Other Information

### DISCLAIMER

For R&D use only. Not for drug, household or other uses.

### WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.



## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278