



Product: Fmoc-Homo-Arg(Pbf)-OH
SDS No: 021120
Revision: 1.0

Safety Data Sheet

Section 1 - Product and Company Information

Product Name: Fmoc-Homo-Arg(Pbf)-OH
Product Code: 021120
Company: CSBio(Shanghai) Ltd
EMERGENCY CONTACTS: Dr. Simon Zhu 0086-21-50781655

Section 2 - Hazards Identification

2.1 Classification of the substance or mixture

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Acute toxicity, oral(Category 4) H302

Skin corrosion/irritation (Category 2) H315

Serious eye damage/eye irritation(Category 2A) H319

Acute toxicity, inhalation (Category 4) H332

Specific target organ toxicity, single exposure;Respiratory system(Category 3) H335

No Resource File

2.2 GHS Label elements, including precautionary statements



Pictogram

Signal word

Warning

Hazard statement(s)

H302 Harmful if swallowed

H315 Causes skin irritation

H319 Causes serious eye irritation

H332 Harmful if inhaled

H335 May cause respiratory irritation

Precautionary statement(s)

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

2.3 Other hazards

Additional precautionary phrases are located throughout the safety data sheet.

Section 3 - Composition/Information on Ingredients

Substance Name: Fmoc-Homo-Arg(Pbf)-OH
CAS No. 1159680-21-3

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Homo-Arg(Pbf)-OH

SDS No: 021120

Revision: 1.0

Formula

$C_{35}H_{42}N_4O_7S$

Section 4 - First Aid Measures

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician

In case of skin contact

Take off contaminated clothing and shoes immediately. Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician. Continue rinsing eyes during transport to hospital.

If swallowed

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5 - Fire Fighting Measures

5.1 Extinguishing media

Suitable extinguishing media: alcohol-resistant foam, dry chemical or carbon dioxide

Unsuitable extinguishing media: no data

5.2 Special hazards arising from the substance or mixture

In combustion toxic fumes may form.

5.3 Advice for firefighters

Wear self-contained breathing apparatus.

Wear protective clothing to prevent contact with skin and eyes.

Section 6 - Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapours, mist or gas.

Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

Do not attempt to take action without suitable protective clothing.

For personal protection see section 8.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains or river. Alert the neighbourhood to the presence of fumes or gas.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Homo-Arg(Pbf)-OH
SDS No: 021120
Revision: 1.0

- 6.3 Methods and materials for containment and cleaning up
Mix with sand or vermiculite. Sweep up and shovel. Transfer to a closable, labelled salvage container for disposal by an appropriate method.
- 6.4 Reference to other sections
For disposal see section 13.

Section 7 - Handling and Storage

- 7.1 Precautions for safe handling
F Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Wash hands thoroughly after handling. Ensure there is sufficient ventilation of the area.
For precautions see section 2.2.
- 7.2 Conditions for safe storage, including any incompatibilities
Keep container tightly closed in a dry and well-ventilated area. Containers which are opened must be carefully resealed and kept upright to prevent leakage. Sealed in dry, $\leq 8^{\circ}\text{C}$
- 7.3 Specific end use(s)
The end use(s) have not been fully determined. The substance is supplied for Research and Development purposes by professionals only.

Section 8 - Exposure Controls / PPE

- 8.1 Control parameters
Components with workplace control parameters
Contains no substances with occupational exposure limit values.
- 8.2 Exposure controls
Appropriate engineering controls
General industrial hygiene practice.
Personal protective equipment
Eye/face protection
Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
Skin protection
Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.
Body Protection
Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
Respiratory protection
Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Homo-Arg(Pbf)-OH
SDS No: 021120
Revision: 1.0

Control of environmental exposure

No special environmental precautions required.

Section 9 - Physical/Chemical Properties

Appearance:	White to off white powder
Molecular Weight:	662.8
pH:	N/A
BP/BP Range:	N/A
MP/MP Range:	N/A
Freezing Point:	N/A
Vapor Pressure:	N/A
Vapor Density:	N/A
Saturated Vapor Conc.	N/A
SG/Density:	N/A
Bulk Density:	N/A
Odor Threshold:	N/A
Volatile%:	N/A
VOC Content:	N/A
Water Content:	≤2.0%
Solvent Content:	N/A
Evaporation Rate:	N/A
Viscosity:	N/A
Surface Tension:	N/A
Partition Coefficient:	N/A
Decomposition Temp.:	N/A
Flash Point:	N/A
Explosion Limits:	N/A
Flammability:	N/A
Autoignition Temp:	N/A
Refractive Index:	N/A
Optical Rotation:	-5.5°±2.0°(C=1, DMF)
Miscellaneous Data:	N/A
Solubility:	N/A

Section 10 - Stability and Reactivity

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Homo-Arg(Pbf)-OH

SDS No: 021120

Revision: 1.0

10.1 Reactivity

No data available

10.2 Chemical stability

Stable under recommended storage conditions.

10.3 Possibility of hazardous reactions

No data available

10.4 Conditions to avoid

No data available

10.5 Incompatible materials

Strong oxidizing agents

10.6 Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NO_x)

Other decomposition products - No data available

In the event of fire: see section 5

Section 11 - Toxicological Information

11.1 Information on toxicological effects

Acute toxicity

LD50 Oral - Rat - male and female - > 5,110 mg/kg

Inhalation: No data available

Dermal: No data available

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Reproductive toxicity

No data available

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Homo-Arg(Pbf)-OH
SDS No: 021120
Revision: 1.0

Specific target organ toxicity - single exposure

No data available

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Additional Information

RTECS: Not available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 - Ecological Information

12.1 Toxicity

No data available

12.2 Persistence and degradability

No data available

12.3 Bioaccumulative potential

No data available

12.4 Mobility in soil

No data available

12.5 Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

12.6 Other adverse effects

No data available

Section 13 - Disposal Considerations

13.1 Waste treatment methods

Product

Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging

Dispose of as unused product.

Section 14 - Transport Information

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

Section 15 - Regulatory Information

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Homo-Arg(Pbf)-OH
SDS No: 021120
Revision: 1.0

SARA 302 Components

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Fmoc-HomoArg(Pbf)-OH CAS-No. 1159680-21-3 Revision Date

New Jersey Right To Know Components

Fmoc-HomoArg(Pbf)-OH CAS-No. 1159680-21-3 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

For R&D use only. Not for drug, household or other uses.

WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278