



Product: Fmoc-Asp-OtBu  
SDS No: 041114  
Revision: 1.0

## Safety Data Sheet

### Section 1 - Product and Company Information

Product Name: Fmoc-Asp-OtBu  
Product Code: 041114  
Company: CSBio(Shanghai) Ltd  
EMERGENCY CONTACTS: Dr. Simon Zhu 0086-21-50781655

### Section 2 - Hazards Identification

#### 2.1 Classification of the substance or mixture

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Acute toxicity, oral(Category 4) H302

Skin corrosion/irritation (Category 2) H315

Serious eye damage/eye irritation(Category 2A ) H319

Specific target organ toxicity, single exposure;Respiratory system(Category 3) H335

No Resource File

#### 2.2 GHS Label elements, including precautionary statements



Pictogram

Signal word

Warning

Hazard statement(s)

H302 Harmful if swallowed

H315 Causes skin irritation

H319 Causes serious eye irritation

H335 May cause respiratory irritation

Precautionary statement(s)

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

#### 2.3 Other hazards

Additional precautionary phrases are located throughout the safety data sheet.

### Section 3 - Composition/Information on Ingredients

Substance Name Fmoc-Asp-OtBu  
CAS No. 129460-09-9  
Formula  $C_{23}H_{25}NO_6$

### Section 4 - First Aid Measures

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Asp-OtBu

SDS No: 041114

Revision: 1.0

#### 4.1 Description of first aid measures

##### General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

##### If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

##### In case of skin contact

Take off contaminated clothing and shoes immediately. Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

##### In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician. Continue rinsing eyes during transport to hospital.

##### If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

#### 4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

#### 4.3 Indication of any immediate medical attention and special treatment needed

No data available

### Section 5 - Fire Fighting Measures

#### 5.1 Extinguishing media

##### Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

#### 5.2 Special hazards arising from the substance or mixture

Carbon oxides, nitrogen oxides (NOx)

#### 5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

#### 5.4 Further information

No data available

### Section 6 - Accidental Release Measures

#### 6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. Do not attempt to take action without suitable protective clothing.

For personal protection see section 8.

#### 6.2 Environmental precautions

No special environmental precautions required.

#### 6.3 Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

#### 6.4 Reference to other sections

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Asp-OtBu  
SDS No: 041114  
Revision: 1.0

For disposal see section 13.

## Section 7 - Handling and Storage

### 7.1 Precautions for safe handling

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.

For precautions see section 2.2.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place ( $\leq 8^{\circ}\text{C}$ ).

### 7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

## Section 8 - Exposure Controls / PPE

### 8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

### 8.2 Exposure controls

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Control of environmental exposure

No special environmental precautions required.

## Section 9 - Physical/Chemical Properties

Appearance: White powder

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Asp-OtBu  
SDS No: 041114  
Revision: 1.0

Molecular Weight:	411.5
pH:	N/A
BP/BP Range:	N/A
MP/MP Range:	N/A
Freezing Point:	N/A
Vapor Pressure:	N/A
Vapor Density:	N/A
Saturated Vapor Conc.	N/A
SG/Density:	N/A
Bulk Density:	N/A
Odor Threshold:	N/A
Volatile%:	N/A
VOC Content:	N/A
Water Content:	≤1.0
Solvent Content:	N/A
Evaporation Rate:	N/A
Viscosity:	N/A
Surface Tension:	N/A
Partition Coefficient:	N/A
Decomposition Temp.:	N/A
Flash Point:	N/A
Explosion Limits:	N/A
Flammability:	N/A
Autoignition Temp:	N/A
Refractive Index:	N/A
Optical Rotation:	-23.0° ± 2.0° (C=1, DMF)
Miscellaneous Data:	N/A
Solubility:	N/A

## Section 10 - Stability and Reactivity

### 10.1 Reactivity

No data available

### 10.2 Chemical stability

Stable under recommended storage conditions.

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Asp-OtBu  
SDS No: 041114  
Revision: 1.0

- 10.3 Possibility of hazardous reactions  
No data available
- 10.4 Conditions to avoid  
No data available
- 10.5 Incompatible materials  
Strong oxidizing agents
- 10.6 Hazardous decomposition products  
Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)  
Other decomposition products - No data available  
In the event of fire: see section 5

## Section 11 - Toxicological Information

- 11.1 Information on toxicological effects
- Acute toxicity  
No data available
- Inhalation: No data available
- Dermal: No data available
- No data available
- Skin corrosion/irritation  
No data available
- Serious eye damage/eye irritation  
No data available
- Respiratory or skin sensitisation  
No data available
- Germ cell mutagenicity  
No data available
- Carcinogenicity  
IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  
ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.  
NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.  
OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
- Reproductive toxicity  
No data available
- Specific target organ toxicity - single exposure  
No data available
- Specific target organ toxicity - repeated exposure  
No data available

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Asp-OtBu  
SDS No: 041114  
Revision: 1.0

Aspiration hazard  
No data available  
Additional Information  
RTECS: Not available

## Section 12 - Ecological Information

- 12.1 Toxicity  
No data available
- 12.2 Persistence and degradability  
No data available
- 12.3 Bioaccumulative potential  
No data available
- 12.4 Mobility in soil  
No data available
- 12.5 Results of PBT and vPvB assessment  
PBT/vPvB assessment not available as chemical safety assessment not required/not conducted
- 12.6 Other adverse effects  
No data available

## Section 13 - Disposal Considerations

- 13.1 Waste treatment methods  
Product  
Offer surplus and non-recyclable solutions to a licensed disposal company.  
Contaminated packaging  
Dispose of as unused product.

## Section 14 - Transport Information

DOT (US)  
Not dangerous goods

IMDG  
Not dangerous goods

IATA  
Not dangerous goods

## Section 15 - Regulatory Information

SARA 302 Components  
No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components  
This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Asp-OtBu  
SDS No: 041114  
Revision: 1.0

No SARA Hazards

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Fmoc-Asp-OtBu CAS-No. 129460-09-9 Revision Date

New Jersey Right To Know Components

Fmoc-Asp-OtBu CAS-No. 129460-09-9 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

## Section 16 - Other Information

### DISCLAIMER

For R&D use only. Not for drug, household or other uses.

### WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278