



Product: H-D-Asp-OH  
SDS No: 041423  
Revision: 1.0

## Safety Data Sheet

### Section 1 - Product and Company Information

Product Name: H-D-Asp-OH  
Product Code: 041423  
Company: CSBio(Shanghai) Ltd  
EMERGENCY CONTACTS: Dr. Simon Zhu 0086-21-50781655

### Section 2 - Hazards Identification

- 2.1 Classification of the substance or mixture  
Not a hazardous substance or mixture.
- 2.2 GHS Label elements, including precautionary statements  
Not a hazardous substance or mixture.
- 2.3 Hazards not otherwise classified (HNOC) or not covered by GHS - none

### Section 3 - Composition/Information on Ingredients

Substance Name H-D-Asp-OH  
CAS No. 1783-96-6  
Formula  $C_4H_7NO_4$

### Section 4 - First Aid Measures

- 4.1 Description of first aid measures  
If inhaled  
If breathed in, move person into fresh air. If not breathing, give artificial respiration.  
In case of skin contact  
Wash off with soap and plenty of water.  
In case of eye contact  
Flush eyes with water as a precaution.  
If swallowed  
Never give anything by mouth to an unconscious person. Rinse mouth with water.
- 4.2 Most important symptoms and effects, both acute and delayed  
The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11
- 4.3 Indication of any immediate medical attention and special treatment needed  
No data available

### Section 5 - Fire Fighting Measures

- 5.1 Extinguishing media  
Suitable extinguishing media  
Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
- 5.2 Special hazards arising from the substance or mixture  
No data available

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: H-D-Asp-OH

SDS No: 041423

Revision: 1.0

### 5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

### 5.4 Further information

No data available

## Section 6 - Accidental Release Measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Avoid dust formation. Avoid breathing vapours, mist or gas.

For personal protection see section 8.

### 6.2 Environmental precautions

No special environmental precautions required.

### 6.3 Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

### 6.4 Reference to other sections

For disposal see section 13.

## Section 7 - Handling and Storage

### 7.1 Precautions for safe handling

Further processing of solid materials may result in the formation of combustible dusts. The potential for combustible dust formation should be taken into consideration before additional processing occurs.

Provide appropriate exhaust ventilation at places where dust is formed.

For precautions see section 2.2.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place ( $\leq 25^{\circ}\text{C}$ ).

### 7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

## Section 8 - Exposure Controls / PPE

### 8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

### 8.2 Exposure controls

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: H-D-Asp-OH  
SDS No: 041423  
Revision: 1.0

use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

#### Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

#### Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

#### Control of environmental exposure

No special environmental precautions required.

### Section 9 - Physical/Chemical Properties

|                        |                          |
|------------------------|--------------------------|
| Appearance:            | White crystalline powder |
| Molecular Weight:      | 133.1                    |
| pH:                    | N/A                      |
| BP/BP Range:           | N/A                      |
| MP/MP Range:           | N/A                      |
| Freezing Point:        | N/A                      |
| Vapor Pressure:        | N/A                      |
| Vapor Density:         | N/A                      |
| Saturated Vapor Conc.  | N/A                      |
| SG/Density:            | N/A                      |
| Bulk Density:          | N/A                      |
| Odor Threshold:        | N/A                      |
| Volatile%:             | N/A                      |
| VOC Content:           | N/A                      |
| Water Content:         | N/A                      |
| Solvent Content:       | N/A                      |
| Evaporation Rate:      | N/A                      |
| Viscosity:             | N/A                      |
| Surface Tension:       | N/A                      |
| Partition Coefficient: | N/A                      |
| Decomposition Temp.:   | N/A                      |
| Flash Point:           | N/A                      |
| Explosion Limits:      | N/A                      |

### CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: H-D-Asp-OH

SDS No: 041423

Revision: 1.0

|                     |                                   |
|---------------------|-----------------------------------|
| Flammability:       | N/A                               |
| Autoignition Temp:  | N/A                               |
| Refractive Index:   | N/A                               |
| Optical Rotation:   | $-25 \pm 2^\circ$ (C=8 in 6N HCl) |
| Miscellaneous Data: | N/A                               |
| Solubility:         | N/A                               |

## Section 10 - Stability and Reactivity

- 10.1 Reactivity
  - No data available
- 10.2 Chemical stability
  - Stable under recommended storage conditions.
- 10.3 Possibility of hazardous reactions
  - No data available
- 10.4 Conditions to avoid
  - No data available
- 10.5 Incompatible materials
  - Strong oxidizing agents
- 10.6 Hazardous decomposition products
  - Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)
  - Other decomposition products - No data available
  - In the event of fire: see section 5

## Section 11 - Toxicological Information

- 11.1 Information on toxicological effects
  - Acute toxicity
    - No data available
  - Inhalation: No data available
  - Dermal: No data available
    - No data available
  - Skin corrosion/irritation
    - No data available
  - Serious eye damage/eye irritation
    - No data available
  - Respiratory or skin sensitisation
    - No data available
  - Germ cell mutagenicity
    - No data available
  - Carcinogenicity
    - IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable,

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: H-D-Asp-OH  
SDS No: 041423  
Revision: 1.0

possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Reproductive toxicity

No data available

Specific target organ toxicity - single exposure

No data available

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Additional Information

RTECS: Not available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

## Section 12 - Ecological Information

### 12.1 Toxicity

No data available

### 12.2 Persistence and degradability

No data available

### 12.3 Bioaccumulative potential

No data available

### 12.4 Mobility in soil

No data available

### 12.5 Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

### 12.6 Other adverse effects

No data available

## Section 13 - Disposal Considerations

### 13.1 Waste treatment methods

Product

Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging

Dispose of as unused product.

## Section 14 - Transport Information

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: H-D-Asp-OH  
SDS No: 041423  
Revision: 1.0

DOT (US)  
Not dangerous goods

IMDG  
Not dangerous goods

IATA  
Not dangerous goods

## Section 15 - Regulatory Information

SARA 302 Components

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

H-D-Asp-OH CAS-No. 1783-96-6 Revision Date

New Jersey Right To Know Components

H-D-Asp-OH CAS-No. 1783-96-6 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

## Section 16 - Other Information

HMIS Rating

Health hazard: 0

Chronic Health Hazard:

Flammability: 0

Physical Hazard 0

NFPA Rating

Health hazard: 0

Fire Hazard: 0

Reactivity Hazard: 0

DISCLAIMER

For R&D use only. Not for drug, household or other uses.

WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a

## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China  
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)  
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA  
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: H-D-Asp-OH

SDS No: 041423

Revision: 1.0

guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.



## CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: [quality@csbiochina.com](mailto:quality@csbiochina.com) • [www.csbiochina.com](http://www.csbiochina.com)

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278