



Product: Fmoc-Cys(Acm)-OH
SDS No: 051101
Revision: 1.0

Safety Data Sheet

Section 1 - Product and Company Information

Product Name: Fmoc-Cys(Acm)-OH
Product Code: 051101
Company: CSBio(Shanghai) Ltd
EMERGENCY CONTACTS: Dr. Simon Zhu 0086-21-50781655

Section 2 - Hazards Identification

2.1 Classification of the substance or mixture

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Acute toxicity, oral(Category 4) H302

Skin corrosion/irritation (Category 2) H315

Serious eye damage/eye irritation(Category 2A) H319

Specific target organ toxicity, single exposure;Respiratory system(Category 3) H335

No Resource File

2.2 GHS Label elements, including precautionary statements



Pictogram

Signal word

Warning

Hazard statement(s)

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P264 Wash skin thoroughly after handling.

P271 Use only outdoors or in a well-ventilated area.

P280 Wear protective gloves/ eye protection/ face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P312 Call a POISON CENTER/doctor if you feel unwell.

P321 Specific treatment (see supplemental first aid instructions on this label).

P332 + P313 If skin irritation occurs: Get medical advice/ attention.

P337 + P313 If eye irritation persists: Get medical advice/ attention.

P362 Take off contaminated clothing and wash before reuse. P

403 + P233 Store in a well-ventilated place. Keep container tightly closed.

P405 Store locked up.

P501 Dispose of contents/ container to an approved waste disposal plant.

2.3 Other hazards

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Cys(Acm)-OH
SDS No: 051101
Revision: 1.0

Additional precautionary phrases are located throughout the safety data sheet.

Section 3 - Composition/Information on Ingredients

Substance Name	Fmoc-Cys(Acm)-OH
CAS No.	86060-81-3
Formula	$C_{21}H_{22}N_2O_5S$

Section 4 - First Aid Measures

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

In case of skin contact

Take off contaminated clothing and shoes immediately. Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician. Continue rinsing eyes during transport to hospital.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5 - Fire Fighting Measures

5.1 Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

5.2 Special hazards arising from the substance or mixture

Carbon oxides, nitrogen oxides (NOx)

5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

5.4 Further information

No data available

Section 6 - Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. Do not attempt to take action without

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Cys(Acm)-OH

SDS No: 051101

Revision: 1.0

suitable protective clothing.

For personal protection see section 8.

6.2 Environmental precautions

No special environmental precautions required.

6.3 Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For disposal see section 13.

Section 7 - Handling and Storage

7.1 Precautions for safe handling

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place ($\leq 8^{\circ}\text{C}$).

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

Section 8 - Exposure Controls / PPE

8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China

Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com

U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA

Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Cys(Acm)-OH
SDS No: 051101
Revision: 1.0

N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Control of environmental exposure

No special environmental precautions required.

Section 9 - Physical/Chemical Properties

Appearance:	White to off white powder
Molecular Weight:	414.5
pH:	N/A
BP/BP Range:	N/A
MP/MP Range:	N/A
Freezing Point:	N/A
Vapor Pressure:	N/A
Vapor Density:	N/A
Saturated Vapor Conc.	N/A
SG/Density:	N/A
Bulk Density:	N/A
Odor Threshold:	N/A
Volatile%:	N/A
VOC Content:	N/A
Water Content:	≤1.0
Solvent Content:	N/A
Evaporation Rate:	N/A
Viscosity:	N/A
Surface Tension:	N/A
Partition Coefficient:	N/A
Decomposition Temp.:	N/A
Flash Point:	N/A
Explosion Limits:	N/A
Flammability:	N/A
Autoignition Temp:	N/A
Refractive Index:	N/A
Optical Rotation:	-32.0° ±3.0° (C=1, MeOH)
Miscellaneous Data:	N/A

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Cys(Acm)-OH
SDS No: 051101
Revision: 1.0

Solubility: N/A

Section 10 - Stability and Reactivity

10.1 Reactivity

No data available

10.2 Chemical stability

Stable under recommended storage conditions.

10.3 Possibility of hazardous reactions

No data available

10.4 Conditions to avoid

No data available

10.5 Incompatible materials

Strong oxidizing agents

10.6 Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)

Other decomposition products - No data available

In the event of fire: see section 5

Section 11 - Toxicological Information

11.1 Information on toxicological effects

Acute toxicity

No data available

Inhalation: No data available

Dermal: No data available

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Cys(Acm)-OH
SDS No: 051101
Revision: 1.0

carcinogen or potential carcinogen by OSHA.
Reproductive toxicity
No data available
Specific target organ toxicity - single exposure
No data available
Specific target organ toxicity - repeated exposure
No data available
Aspiration hazard
No data available
Additional Information
RTECS: Not available

Section 12 - Ecological Information

- 12.1 Toxicity
No data available
- 12.2 Persistence and degradability
No data available
- 12.3 Bioaccumulative potential
No data available
- 12.4 Mobility in soil
No data available
- 12.5 Results of PBT and vPvB assessment
PBT/vPvB assessment not available as chemical safety assessment not required/not conducted
- 12.6 Other adverse effects
No data available

Section 13 - Disposal Considerations

- 13.1 Waste treatment methods
Product
Offer surplus and non-recyclable solutions to a licensed disposal company.
Contaminated packaging
Dispose of as unused product.

Section 14 - Transport Information

DOT (US)
Not dangerous goods

IMDG
Not dangerous goods

IATA
Not dangerous goods

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Fmoc-Cys(Acm)-OH
SDS No: 051101
Revision: 1.0

Section 15 - Regulatory Information

SARA 302 Components

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Fmoc-Cys(Acm)-OH CAS-No. 86060-81-3 Revision Date

New Jersey Right To Know Components

Fmoc-Cys(Acm)-OH CAS-No. 86060-81-3 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

DISCLAIMER

For R&D use only. Not for drug, household or other uses.

WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278