



## SPECIFICATION

| Tests                               | Specifications                    |
|-------------------------------------|-----------------------------------|
| Appearance                          | White to off white powder         |
| Purity (HPLC) (%)                   | ≥98.0                             |
| Specific Rotation $[\alpha]_D^{20}$ | +9.0°±2.0°(C=1 in DMF)            |
| Water Content(K.F.) (%)             | ≤1.0                              |
| Clarity of solution                 | 5 mmol in 10ml DMF clear solution |
| Identification(MS)                  | In accordance with the standard   |
| Identification(NMR)                 | In accordance with the standard   |

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information