



## SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Assay(Titration) (%)	98.0-102.0
Specific Rotation $[\alpha]_D^{20}$	$+7.0^{\circ} \pm 1.0^{\circ}$ (C=2, MeOH)
Loss on drying(%)	$\leq 1.0$
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information