



SPECIFICATION

Tests	Specifications
Appearance	White crystalline powder
Assay(%)	98.0-102.0
Transmittance	≥ 95
Chloride (Cl) (%)	≤ 0.04
Sulphate (SO ₄) (%)	≤ 0.05
Ammonium(NH ₄) (%)	≤ 0.02
Iron(Fe)	≤ 30.0 ppm
Heavy metals (Pb)	≤ 10.0 ppm
Arsenic(As ₂ O ₃)	≤ 2.0 ppm
Loss on drying(%)	≤ 0.5
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information