



SPECIFICATION

Tests	Specifications
Appearance	white powder
Purity(UPLC) (%)	≥99.0
Optical Purity(%)	≤0.5 L-enantiomer
Specific Rotation $[\alpha]_D^{20}$	+4.5°±1.5°(C=1, DMF)
Clarity of solution	0.3 gram in 2ml DMF clear solution
Water Content(K.F.) (%)	≤2.0
Identification (MS)	In accordance with the standard
Identification(NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information