



SPECIFICATION

Tests	Specifications
Appearance	White powder
Purity(HPLC) (%)	≥98.0
Specific Rotation $[\alpha]_D^{20}$	-3.5°±1.0°(C=2, EtOH)
Clarity of solution	0.3 gram in 2ml EtOH clear solution
Melting point	113.0-123.0°C
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information