



SPECIFICATION

Tests	Specifications
Appearance	White powder
Assay (%)	99.0-102.0
Loss on drying (%)	11.5-12.5
Sulphate(SO ₄) (%)	≤0.02
Ammonium(NH ₄) (%)	≤0.1
Chloride(Cl) (%)	≤0.02
Heavy Metals(Pb)	≤10ppm
Arsenic(As ₂ O ₃)	≤2ppm
Iron(Fe)	≤10ppm
Residue on ignition (%)	≤0.1
Identification(MS)	In accordance with the standard
Identification(NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information