



## SPECIFICATION

Tests	Specifications
Appearance	White powder
Purity(HPLC) (%)	≥98.0
Specific Rotation $[\alpha]_D^{20}$	-27.0° ± 2.0° (C=1, DMF)
Melting Point	110.0-115.0°C
Clarity of solution	0.3 gram in 2ml DMF clear solution
Identification(MS)	In accordance with the standard
Identification(NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information