



## SPECIFICATION

| Tests                               | Specifications                  |
|-------------------------------------|---------------------------------|
| Appearance                          | White to off white powder       |
| Purity(HPLC) (%)                    | ≥95.0                           |
| Specific Rotation $[\alpha]_D^{20}$ | -24.5° ± 2.5° (C=1, DMF)        |
| Identification(MS)                  | In accordance with the standard |
| Identification(NMR)                 | In accordance with the standard |

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information