



SPECIFICATION

Tests	Specifications
Appearance	White powder
Purity(HPLC) (%)	≥98.0
Optical Purity (%)	≤0.5 D-enantiomer
Specific Rotation $[\alpha]_D^{20}$	-14.0°±2.0° (C=1, MeOH)
Melting Point	80.0-90.0°C
Clarity of solution	0.3 gram in 2ml MeOH clear solution
Water Content(K.F.)	≤1.0
Identification (MS) (%)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information