



SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Assay (%)	98.0-102.0
Specific Rotation $[\alpha]_D^{20}$	$+9.0^{\circ} \pm 2.0^{\circ}$ (C =2, H ₂ O)
Melting Point	240.0-250.0°C
Water Content(K.F.) (%)	≤1.0
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information