



=

SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Purity(HPLC) (%)	≥98.0
Specific Rotation $[\alpha]_D^{20}$	+25.0°±3.0°(C=1, MeOH)
Water Content(K.F.)	≤2.0
Clarity of solution	0.3 gram in 2ml DMF clear solution
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information