



## SPECIFICATION

| Tests                               | Specifications  |
|-------------------------------------|---|
| Appearance                          | White crystalline powder                                |
| Purity(HPLC) (%)                    | ≥97.0   |
| Specific Rotation $[\alpha]_D^{20}$ | -23.0° ±3.0° (C = 1 , EtOH)                             |
| Clarity of solution                 | 0.5g in 10 ml MeOH clear to very slightly hazy solution |
| Identification (MS)                 | In accordance with the standard                         |
| Identification (NMR)                | In accordance with the standard                         |

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information