



SPECIFICATION

Tests	Specifications
Appearance	White crystalline powder
Purity(HPLC) (%)	≥98.5
Impurity(%)	Fmoc-β-Ala-OH≤0.5
Melting Point	182.0-188.0°C
Loss on drying(%)	≤0.5
Clarity of solution	5 mmol in 10ml DMF clear solution
Water Content(K.F.) (%)	≤3.0
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information