



SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Purity(HPLC) (%)	≥98.0
Specific Rotation $[\alpha]_D^{20}$	-2.0° ±1.0° (C =1 , MeOH)
Clarity of solution	5 mmol in 10ml DMF clear solution
Water Content(K.F.) (%)	≤2.0
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information