



SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Purity(HPLC) (%)	≥99.0
Specific Rotation $[\alpha]_D^{20}$	+29.0° ±2.0° (C =1, DMF)
Optical Purity(%)	≤0.3 L-enantiomer
Clarity of solution	0.3 gram in 2ml DMF clear solution
Water Content(K.F.) (%)	≤2.5
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information