



SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Assay (%)	98.5-100.5
Chloride (Cl) (%)	≤ 0.02
Ammonium(NH ₄) (%)	≤ 0.02
Iron(Fe)	$\leq 30\text{ppm}$
Heavy metals (Pb)	$\leq 10\text{ppm}$
Arsenic(As ₂ O ₃)	$\leq 1\text{ppm}$
Loss on drying (%)	≤ 0.2
Residue on ignition (%)	≤ 0.15
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information