



SPECIFICATION

Tests	Specifications
Appearance	White powder
Purity(HPLC) (%)	≥98.0
Specific Rotation $[\alpha]_D^{20}$	-65.0°±2.5° (C =1, CHCl ₃)
Melting point	92.0-96.0°C
Clarity of solution	0.3 gram in 2ml DMF clear solution
Water Content(K.F.) (%)	≤1.0
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information