



## SPECIFICATION

| Tests                | Specifications                  |
|----------------------|---------------------------------|
| Appearance           | White crystalline powder        |
| Purity(HPLC) (%)     | ≥99.0                           |
| Melting Point        | 196.0-205.0°C                   |
| Loss on drying (%)   | ≤0.5                            |
| Clarity of solution  | 1 gram in 10ml ACN clear        |
| Identification (MS)  | In accordance with the standard |
| Identification (NMR) | In accordance with the standard |

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information