



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

Safety Data Sheet

Section 1 - Product and Company Information

Product Name: Boc Anhydride
Product Code: 321004
Company: CSBio(Shanghai) Ltd
EMERGENCY CONTACTS: Dr. Simon Zhu +86-21-50781655

Section 2 - Hazards Identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008

Flammable liquids (Category 3), H226

Acute toxicity, Inhalation (Category 1), H330

Skin irritation (Category 2), H315 Serious eye damage (Category 1), H318

Skin sensitisation (Category 1), H317

Specific target organ toxicity - single exposure (Category 3), Respiratory system, H335

For the full text of the H-Statements mentioned in this Section, see Section 16.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram



Signal word

Danger

Hazard statement(s)

H226 Flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H330 Fatal if inhaled.

H335 May cause respiratory irritation.

Precautionary statement(s)

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P260 Do not breathe dust/ fume/ gas/ mist/ vapours/ spray.

P280 Wear eye protection/ face protection.

P304 + P340 + P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor.

P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. P370 + P378 In case of fire: Use dry powder or dry sand to extinguish.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

Supplemental Hazard Statements none

2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Section 3 - Composition/Information on Ingredients

Substance Name	Boc Anhydride
CAS No.	24424-99-5
Formula	C ₁₀ H ₁₈ O ₅

Section 4 - First Aid Measures

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5 - Fire Fighting Measures

5.1 Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

5.2 Special hazards arising from the substance or mixture

Carbon oxides

5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

5.4 Further information

Use water spray to cool unopened containers.

Section 6 - Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Remove all

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

sources of ignition. Evacuate personnel to safe areas. Beware of vapours accumulating to form explosive concentrations. Vapours can accumulate in low areas.

For personal protection see section 8.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

6.3 Methods and materials for containment and cleaning up

Contain spillage, and then collect with an electrically protected vacuum cleaner or by wetbrushing and place in container for disposal according to local regulations (see section 13).

6.4 Reference to other sections

For disposal see section 13.

Section 7 - Handling and Storage

7.1 Precautions for safe handling

Avoid contact with skin and eyes. Avoid inhalation of vapour or mist. Keep away from sources of ignition - No smoking. Take measures to prevent the build up of electrostatic charge.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Store in cool place. Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Recommended storage temperature $\leq 8^{\circ}\text{C}$

Avoid heating above: 40°C Store under inert gas. Moisture sensitive. Handle and open container with care.

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

Section 8 - Exposure Controls / PPE

8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Control of environmental exposure

No special environmental precautions required.

Section 9 - Physical/Chemical Properties

Appearance:	White crystals or Colorless liquid
Molecular Weight:	218.3
pH:	N/A
BP/BP Range:	N/A
MP/MP Range:	N/A
Freezing Point:	N/A
Vapor Pressure:	N/A
Vapor Density:	N/A
Saturated Vapor Conc.	N/A
SG/Density:	N/A
Bulk Density:	N/A
Odor Threshold:	N/A
Volatile%:	N/A
VOC Content:	N/A
Water Content:	≤0.2%
Solvent Content:	N/A
Evaporation Rate:	N/A
Viscosity:	N/A
Surface Tension:	N/A
Partition Coefficient:	N/A
Decomposition Temp.:	N/A
Flash Point:	N/A
Explosion Limits:	N/A
Flammability:	N/A
Autoignition Temp:	N/A

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

Refractive Index: N/A
Optical Rotation: N/A
Miscellaneous Data: N/A
Solubility: N/A

Section 10 - Stability and Reactivity

- 10.1 Reactivity
No data available
- 10.2 Chemical stability
Stable under recommended storage conditions.
- 10.3 Possibility of hazardous reactions
No data available
- 10.4 Conditions to avoid
Heat, flames and sparks.
- 10.5 Incompatible materials
Oxidizing agents, Reducing agents, Bases, Avoid contact with:, Acids
- 10.6 Hazardous decomposition products
Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NO_x)
Other decomposition products - No data available
In the event of fire: see section 5

Section 11 - Toxicological Information

- 11.1 Information on toxicological effects
Acute toxicity
LD50 Oral - Rat - > 5.000 mg/kg
Remarks: (External MSDS)
LC50 Inhalation - Rat - 4 h - 0,1 mg/l
Remarks: (RTECS) LD50 Dermal - Rabbit - > 2.000 mg/kg
Remarks: (External MSDS)
Skin corrosion/irritation
Serious eye damage/eye irritation
Respiratory or skin sensitisation
Germ cell mutagenicity
No data available
Ames test
Salmonella typhimurium
Result: negative
(National Toxicology Program)
Carcinogenicity
IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

Reproductive toxicity

No data available

Specific target organ toxicity - single exposure

No data available

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Additional Information

RTECS: HT0230000

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 - Ecological Information

12.1 Toxicity

No data available

12.2 Persistence and degradability

No data available

12.3 Bioaccumulative potential

No data available

12.4 Mobility in soil

No data available

12.5 Results of PBT and vPvB assessment

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

12.6 Other adverse effects

No data available

Section 13 - Disposal Considerations

13.1 Waste treatment methods

Product

Burn in a chemical incinerator equipped with an afterburner and scrubber but exert extra care in igniting as this material is highly flammable. Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging

Dispose of as unused product.

Section 14 - Transport Information

14.1 UN number

ADR/RID: 2929 IMDG: 2929 IATA: 2929

14.2 UN proper shipping name

ADR/RID: TOXIC LIQUID, FLAMMABLE, ORGANIC, N.O.S. (Di(tert-butyl) carbonate)

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: Boc Anhydride
SDS No: 321004
Revision: 1.0

IMDG: TOXIC LIQUID, FLAMMABLE, ORGANIC, N.O.S. (Di(tert-butyl) carbonate)

IATA: Toxic liquid, flammable, organic, n.o.s. (Di(tert-butyl) carbonate)

14.3 Transport hazard class(es)

ADR/RID: 6.1 (3) IMDG: 6.1 (3) IATA: 6.1 (3)

14.4 Packaging group

ADR/RID: II IMDG: II IATA: II

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

No data available

Section 15 - Regulatory Information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This safety datasheet complies with the requirements of Regulation (EC) No. 1907/2006.

15.2 Chemical safety assessment

For this product a chemical safety assessment was not carried out

Section 16 - Other Information

Full text of H-Statements referred to under sections 2 and 3.

H226 Flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H319 Causes serious eye irritation.

H330 Fatal if inhaled.

H335 May cause respiratory irritation.

DISCLAIMER

For R&D use only. Not for drug, household or other uses.

WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278