



SPECIFICATION

Tests	Specifications
Appearance	White crystalline powder
Assay (%)	98.5-101.5
Specific Rotation $[\alpha]_D^{20}$	$-11.0^\circ \pm 2.0^\circ$ (C =5, 1N HCl)
Heavy metals (Pb)	$\leq 15\text{ppm}$
Loss on drying (%)	≤ 0.5
Residue on ignition (%)	≤ 0.4

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information