



## SPECIFICATION

Tests	Specifications
Appearance	White crystals or crystalline powder
Assay (%)	98.5-101.5
Specific Rotation $[\alpha]_D^{20}$	$+15.0^{\circ} \pm 2.0^{\circ}$ (C=4, 6N HCl)
Heavy metals(Pb)	$\leq 15\text{ppm}$
Loss on drying (%)	$\leq 0.5$
Residue on ignition (%)	$\leq 0.4$

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information