



SPECIFICATION

Tests	Specifications
Appearance	White crystals to crystalline powder
Assay (%)	98.5-101.5
Specific Rotation $[\alpha]_D^{20}$	$+12.0^{\circ} \pm 3.0^{\circ}$ (C = 1, 6N HCl)
Heavy metals (Pb)	$\leq 15\text{ppm}$
Loss on drying (%)	≤ 0.5
Residue on ignition (%)	≤ 0.3

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information