



SPECIFICATION

Tests	Specifications
Appearance	White crystalline or crystalline powder
Assay (%)	98.5-101.5
Specific Rotation $[\alpha]_D^{20}$	$+24.0^{\circ} \pm 2.0^{\circ}$ (C=2, 6N HCl)
Heavy metals(by Pb)	≤ 10 ppm
Loss on drying (%)	≤ 0.5
Residue on ignition (%)	≤ 0.1

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information