



SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Purity(UPLC) (%)	≥95.0
Optical Purity(%)	≤0.5 L-enantiomer
Specific Rotation $[\alpha]_D^{20}$	+14.5°±2.0°(C=1, DMF)
Clarity of solution	5mmol in 10ml DMF clear solution
Water Content(K.F.) (%)	≤2.0
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information