



SPECIFICATION

Tests	Specifications
Appearance	White to off white powder
Purity(HPLC) (%)	≥90.0
Specific Rotation $[\alpha]_D^{20}$	-3.0°±2.0°(C=4, DMF)
Water Content (K.F.) (%)	≤2.0
Clarity of solution	0.3 gram in 2ml DMF clear solution
Identification (MS)	In accordance with the standard
Identification (NMR)	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information