



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

Safety Data Sheet

Section 1 - Product and Company Information

Product Name: CI-HOBt
Product Code: 311004
Company: CSBio(Shanghai) Ltd
EMERGENCY CONTACTS: Dr. Simon Zhu 0086-21-50781655

Section 2 - Hazards Identification

- 2.1 Classification of the substance or mixture
Acute toxicity - Category 4, Oral
- 2.2 GHS Label elements, including precautionary statements



Pictogram

Signal word

Warning

Hazard statement(s)

H203 Explosive; fire, blast or projection hazard

H302 Harmful if swallowed

Precautionary statement(s)

Prevention

P264 Wash ... thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

Response

P301+P317 IF SWALLOWED: Get medical help.

P330 Rinse mouth.

Storage

None

Disposal

P501 Dispose of contents/container to an appropriate treatment and disposal facility in accordance with applicable laws and regulations, and product characteristics at time of disposal.

- 2.3 Other hazards

Additional precautionary phrases are located throughout the safety data sheet.

Section 3 - Composition/Information on Ingredients

Substance Name: CI-HOBt
CAS No.: 26198-19-6
Formula: $C_6H_4ClN_3O$

Section 4 - First Aid Measures

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

In case of skin contact

Take off contaminated clothing and shoes immediately. Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician. Continue rinsing eyes during transport to hospital.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5 - Fire Fighting Measures

5.1 Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

5.2 Special hazards arising from the substance or mixture

Carbon oxides, nitrogen oxides (NOx)

5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

5.4 Further information

No data available

Section 6 - Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. Do not attempt to take action without suitable protective clothing.

For personal protection see section 8.

6.2 Environmental precautions

No special environmental precautions required.

6.3 Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

For disposal see section 13.

Section 7 - Handling and Storage

7.1 Precautions for safe handling

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Keep in a dry place ($\leq 25^{\circ}\text{C}$, away from light).

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

Section 8 - Exposure Controls / PPE

8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Appropriate engineering controls

General industrial hygiene practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Control of environmental exposure

No special environmental precautions required.

Section 9 - Physical/Chemical Properties

Appearance: White to off white crystalline powder

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

Molecular Weight:	169.6
pH:	N/A
BP/BP Range:	N/A
MP/MP Range:	190.0-200.0°C
Freezing Point:	N/A
Vapor Pressure:	N/A
Vapor Density:	N/A
Saturated Vapor Conc.	N/A
SG/Density:	N/A
Bulk Density:	N/A
Odor Threshold:	N/A
Volatile%:	N/A
VOC Content:	N/A
Water Content:	N/A
Solvent Content:	N/A
Evaporation Rate:	N/A
Viscosity:	N/A
Surface Tension:	N/A
Partition Coefficient:	N/A
Decomposition Temp.:	N/A
Flash Point:	N/A
Explosion Limits:	N/A
Flammability:	N/A
Autoignition Temp:	N/A
Refractive Index:	N/A
Optical Rotation:	N/A
Miscellaneous Data:	N/A
Solubility:	N/A

Section 10 - Stability and Reactivity

10.1 Reactivity

No data available

10.2 Chemical stability

Stable under recommended storage conditions.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

- 10.3 Possibility of hazardous reactions
No data available
- 10.4 Conditions to avoid
No data available
- 10.5 Incompatible materials
Strong oxidizing agents
- 10.6 Hazardous decomposition products
Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)
Other decomposition products - No data available
In the event of fire: see section 5

Section 11 - Toxicological Information

- 11.1 Information on toxicological effects
- Acute toxicity
No data available
- Inhalation: No data available
- Dermal: No data available
- No data available
- Skin corrosion/irritation
No data available
- Serious eye damage/eye irritation
No data available
- Respiratory or skin sensitisation
No data available
- Germ cell mutagenicity
No data available
- Carcinogenicity
IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.
NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
- Reproductive toxicity
No data available
- Specific target organ toxicity - single exposure
No data available
- Specific target organ toxicity - repeated exposure
No data available

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

Aspiration hazard
No data available
Additional Information
RTECS: Not available

Section 12 - Ecological Information

- 12.1 Toxicity
No data available
- 12.2 Persistence and degradability
No data available
- 12.3 Bioaccumulative potential
No data available
- 12.4 Mobility in soil
No data available
- 12.5 Results of PBT and vPvB assessment
PBT/vPvB assessment not available as chemical safety assessment not required/not conducted
- 12.6 Other adverse effects
No data available

Section 13 - Disposal Considerations

- 13.1 Waste treatment methods
Product
Offer surplus and non-recyclable solutions to a licensed disposal company.
Contaminated packaging
Dispose of as unused product.

Section 14 - Transport Information

DOT (US)
Not dangerous goods

IMDG
Not dangerous goods

IATA
Not dangerous goods

Section 15 - Regulatory Information

SARA 302 Components
No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components
This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278



Product: CI-HOBt
SDS No: 311004
Revision: 1.0

No SARA Hazards

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

CI-HOBt CAS-No. 26198-19-6 Revision Date

New Jersey Right To Know Components

CI-HOBt CAS-No. 26198-19-6 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

DISCLAIMER

For R&D use only. Not for drug, household or other uses.

WARRANTY

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CSBio (Shanghai) Ltd., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.

CSBio Co.

Pilot Free Trade Zone, 85 Faladi Road, Building 2, Room 502, Shanghai 201210, China
Tel: +86-21-50781655 • Fax: +86-21-50781652 • Email: quality@csbiochina.com • www.csbiochina.com
U.S. Headquarter: 20 Kelly Court, Menlo Park, CA 94025 USA
Tel: +1-650-3221111 • Fax: +1-650-3222278