



#### SPECIFICATION

Tests	Specifications
Appearance	White to off white crystalline powder
Assay (%)	98.5-101.0
Specific Rotation $[\alpha]_D^{20}$	$-27.0^{\circ} \pm 2.0^{\circ}$ (C=8, 6N HCl)
Heavy metals	$\leq 10.0\text{ppm}$
Loss on drying (%)	$\leq 0.5$
Residue on ignition (%)	$\leq 0.2$

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information