



Product Specification

ORIGINAL

Product Information

Product Name:	H-Asp(OtBu)-2-Chlorotrityl Resin	Catalog No.:	262007
Molecular Weight:	N/A	CAS No.:	N/A
Shelf-life	24months	Molecular Formula:	N/A
Recommended Long Term Storage Condition:	Cool dry place ($\leq 8^{\circ}\text{C}$)		
Grade:	Standard Grade		

SPECIFICATION

Tests	Specifications
Appearance	Yellow beads
Cross-Linking	1%DVB
Bead Size	100-200mesh
Loading	0.3-0.8mmol/g
IR Spectrum	In accordance with the standard

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information

