



## Product Specification

ORIGINAL

### Product Information

|                                          |                                             |                    |        |
|------------------------------------------|---------------------------------------------|--------------------|--------|
| Product Name:                            | H-Glu(OtBu)-2-Chlorotrityl Resin            | Catalog No.:       | 262016 |
| Molecular Weight:                        | N/A                                         | CAS No.:           | N/A    |
| Shelf-life                               | 24months                                    | Molecular Formula: | N/A    |
| Recommended Long Term Storage Condition: | Cool dry place ( $\leq 8^{\circ}\text{C}$ ) |                    |        |
| Grade:                                   | Standard Grade                              |                    |        |

### SPECIFICATION

| Tests         | Specifications                  |
|---------------|---------------------------------|
| Appearance    | Yellow beads                    |
| Cross-Linking | 1%DVB                           |
| Bead Size     | 100-200mesh                     |
| Loading       | 0.3-0.8mmol/g                   |
| IR Spectrum   | In accordance with the standard |

Note: CS Bio reserves the right to revise the above information without any written consent to any users on any intent use of this information

